

Summary

Name	ANGPTL7/Angiopietin-related Protein 7
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Angiopietin-related Protein 7 is produced by our Mammalian expression system and the target gene encoding Gln27-Pro346 is expressed with a 10His tag at the C-terminus.
Accession #	O43827
Host	Human Cells
Species	Human
Predicted Molecular Mass	38.4 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 6% Trehalose, 4% Mannitol, 100mM NaCl, 0.05% Tween 80, pH6.0.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human ANGPTL7 (C-10His)
Catalog #: PHH0065



Alternative Names

Angiopoietin-related protein 7; Angiopoietin-like factor; Angiopoietin-like protein 7; Cornea-derived transcript 6 protein; ANGPTL7; CDT6

Background

Angiopoietin-like 7 (ANGPTL7) is a secreted glycoprotein that is structurally related to the angiopoietins. Members of this protein family contain an N-terminal coiled coil domain and a C-terminal fibrinogen-like domain. ANGPTL7 shares 89% aa sequence identity with mouse and rat ANGPTL7. It is secreted as a 45-50kDa monomer that forms disulfide-linked homotrimers and tetramers via the coiled coil domain. ANGPTL7 is expressed in the corneal stroma, trabecular meshwork, and sclera and is elevated in glaucoma aqueous humor. Its production is up-regulated in trabecular meshwork cells by glucocorticoids and TGF- β and in cartilage by TNF- α . Overexpression of ANGPTL7 in trabecular meshwork cells inhibits the production of collagen and proteoglycans. When overexpressed in tumor cells it promotes collagen and proteoglycan deposition but inhibits tumor xenograft progression and tumor angiogenesis.

Note

For Research Use Only , Not for Diagnostic Use.