

Product Name: Recombinant Human BCAN (C-6His)
Catalog #: PHH0169




Summary

Name	Brevican core protein/BCAN
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Brevican Core Protein is produced by our Mammalian expression system and the target gene encoding Asp23-Pro911 is expressed with a 6His tag at the C-terminus.
Accession #	AAH09117.1
Host	Human Cells
Species	Human
Predicted Molecular Mass	97.75 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human BCAN (C-6His)
Catalog #: PHH0169



Alternative Names

Brevican Core Protein; Brain-Enriched Hyaluronan-Binding Protein; BEHAB; Chondroitin Sulfate Proteoglycan 7; BCAN; BEHAB; CSPG7

Background

Brevican Core Protein (BCAN) is a secreted protein that belongs to the aggrecan/versican proteoglycan family. BCAN contains one C-type lectin domain, one EGF-like domain, one Ig-like V-type domain, one Sushi (CCP/SCR) domain and two Link domains. BCAN may play a role in the terminally differentiating and the adult nervous system during postnatal development. BCAN could stabilize interactions between hyaluronan (HA) and brain proteoglycans.

Note

For Research Use Only , Not for Diagnostic Use.