


**Product Name: Recombinant Human CD27L (C-Fc)**  
**Catalog #: PHH0368**



## Summary

<b>Name</b>	CD27L/CD70/CD27 Ligand/TNFSF7
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	<1 EU/μg as determined by LAL test.
<b>Construction</b>	Recombinant Human CD27 Ligand/CD70 is produced by our Mammalian expression system and the target gene encoding Gln45-Pro193 is expressed with a human IgG1 Fc tag at the C-terminus.
<b>Accession #</b>	P32970
<b>Host</b>	Human Cells
<b>Species</b>	Human
<b>Predicted Molecular Mass</b>	43.6 KDa
<b>Formulation</b>	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Reconstitution</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE image



## Background

<b>Alternative Names</b>	CD70 antigen; CD70; CD27 ligand; Tumor necrosis factor ligand superfamily
--------------------------	---

**Product Name: Recombinant Human CD27L (C-Fc)**  
**Catalog #: PHH0368**



---

**Background**

member 7; CD27LG; TNFSF7; CD27L

CD27 Ligand (CD70) is a 30 kDa type 2 transmembrane glycoprotein that belongs to the tumor necrosis factor (TNF) ligand family. Human CD27 Ligand extracellular region is 64% aa identical to mouse and rat CD27 Ligand extracellular regions. This cytokine is a ligand for TNFRSF27/CD27. It is a surface antigen on activated, but not on resting, T and B lymphocytes. It induces proliferation of costimulated T cells, enhances the generation of cytolytic T cells, and contributes to T cell activation. This cytokine is also reported to play a role in regulating B-cell activation, cytotoxic function of natural killer cells, and immunoglobulin synthesis.

**Note**

For Research Use Only , Not for Diagnostic Use.