

Product Name: Recombinant Human CD99L2 (C-6His)
Catalog #: PHH0383




Summary

Name	CD99 antigen-like protein 2/CD99L2
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human CD99 Antigen-Like Protein 2 is produced by our Mammalian expression system and the target gene encoding Asp26-Ala188 is expressed with a 6His tag at the C-terminus.
Accession #	Q8TCZ2
Host	Human Cells
Species	Human
Predicted Molecular Mass	18.4 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human CD99L2 (C-6His)
Catalog #: PHH0383



Alternative Names

CD99 Antigen-Like Protein 2; MIC2-Like Protein 1; CD99; CD99L2; MIC2L1

Background

CD99 Antigen-Like Protein 2 (CD99L2) belongs to the CD99 family. CD99L2 is a single-pass type I membrane protein and expressed in many tissues, with low expression in thymus. CD99L2 plays a role in a late step of leukocyte extravasation helping cells to overcome the endothelial basement membrane. CD99L2 and CD99 are involved in trans-endothelial migration of neutrophils in vitro and in the recruitment of neutrophils into inflamed peritoneum. A similar protein in mouse functions as an adhesion molecule during leukocyte extravasation. Alternate splicing results in multiple transcript variants.

Note

For Research Use Only , Not for Diagnostic Use.