

Product Name: Recombinant Human CYB5B (C-6His)
Catalog #: PHH0505




Summary

Name	Cytochrome b5 type B/CYB5B
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Cytochrome B5 Type B is produced by our Mammalian expression system and the target gene encoding Lys12-Cys118 is expressed with a 6His tag at the C-terminus.
Accession #	AAH04373.1
Host	Human Cells
Species	Human
Predicted Molecular Mass	13.1 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human CYB5B (C-6His)
Catalog #: PHH0505



Alternative Names

Cytochrome b5 outer mitochondrial membrane isoform; OMB5

Background

Cytochrome b5 type B (CYB5B) is a membrane of the cytochrome b5 family. It contains 1 cytochrome b5 heme-binding domain. Cytochrome b5 is a membrane bound hemoprotein which function as an electron carrier for several membrane bound oxygenases. In the mitochondrion of eukaryotes and in aerobic prokaryotes, cytochrome b is a component of respiratory chain complex III also known as the bc1 complex or ubiquinol-cytochrome c reductase.

Note

For Research Use Only , Not for Diagnostic Use.