

Product Name: Recombinant Human EFNA3 (C-6His)
Catalog #: PHH0588




Summary

Name	Ephrin-A3/EFNA3
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Ephrin-A3 is produced by our Mammalian expression system and the target gene encoding Gln23-Ser211 is expressed with a 6His tag at the C-terminus.
Accession #	P52797
Host	Human Cells
Species	Human
Predicted Molecular Mass	22.25 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human EFNA3 (C-6His)
Catalog #: PHH0588



Alternative Names

Ephrin-A3; EFL-2; EHK1 Ligand; EHK1-L; EPH-Related Receptor Tyrosine Kinase Ligand 3; LERK-3; EFNA3; EFL2; EPLG3; LERK3

Background

Ephrins-A3 belongs the Ephrins ligand family which involved in a variety of biological processes, especially in the nervous system and in erythropoiesis. It is shown that Ephrin-A3 is expressed in brain, skeletal muscle, spleen, thymus, prostate, testis, ovary, small intestine, and peripheral blood leukocytes. Ephrin-A3 has a GPI anchor following the extracellular sequence and a signal sequence of 22 amino acids. Ephrin-A3 can bind EphA2, EphA3, EphA4, EphA5, EphA6, EphA7, EphA8 and EphB1. Furthermore, it is associated with tumor growth and metastasis.

Note

For Research Use Only , Not for Diagnostic Use.