

Product Name: Recombinant Human ERMAP (C-6His)
Catalog #: PHH0604




Summary

Name	Erythroid membrane-associated protein/ERMAP
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Erythroid Membrane-Associated Protein is produced by our Mammalian expression system and the target gene encoding His30-Ala155 is expressed with a 6His tag at the C-terminus.
Accession #	Q96PL5
Host	Human Cells
Species	Human
Predicted Molecular Mass	14.79 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human ERMAP (C-6His)
Catalog #: PHH0604



Alternative Names

Erythroid Membrane-Associated Protein; hERMAP; Radin Blood Group Antigen; Scianna Blood Group Antigen; ERMAP; RD; SC

Background

Human Erythroid Membrane-Associated Protein (ERMAP) is a cell surface transmembrane protein that belongs to the immunoglobulin superfamily. It is highly expressed in bone marrow and to a lower extent in leukocytes, thymus, lymph node and spleen. ERMAP contains 1 B30.2/SPRY domain and 1 Ig-like V-type (immunoglobulin-like) domain. It may serve as an erythroid cell receptor, possibly as a mediator of cell adhesion. ERMAP is responsible for the Scianna/Radin blood group system. Two transcript variants encoding the same protein have been found for this gene ERMAP.

Note

For Research Use Only , Not for Diagnostic Use.