

Product Name: Recombinant Human INSL4 (C-6His)
Catalog #: PHH0957




Summary

Name	INSL4/Placentin
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Early Placenta Insulin-Like Peptide is produced by our Mammalian expression system and the target gene encoding Ala26-Thr139 is expressed with a 6His tag at the C-terminus.
Accession #	Q14641
Host	Human Cells
Species	Human
Predicted Molecular Mass	13.6 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human INSL4 (C-6His)
Catalog #: PHH0957



Alternative Names

Early Placenta Insulin-Like Peptide; EPIL; Insulin-Like Peptide 4; Placentin; INSL4

Background

Early Placenta Insulin-Like Peptide (INSL4) belongs to the insulin family. INSL4 is expressed in the early placental cytotrophoblast and syncytiotrophoblast. INSL4 is a secreted protein and a precursor that undergoes post-translational cleavage to produce 3 polypeptide chains, A-C, that form tertiary structures composed of either all three chains, or just the A and B chains. INSL4 plays an important role in the development of trophoblast and in the regulation of bone formation.

Note

For Research Use Only , Not for Diagnostic Use.