



Product Name: Recombinant Human F8A3,F8A1,F8A2 Protein
Catalog #: PHH10052



Summary

Name	HAP40_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P23610
Host	HEK293
Species	Human
Molecular Mass	66 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DFLARYRLVS NKLKKRFLRK PNVAEAGEQF GQLGRELRAQ ECLPYAAWCQ LAVARCQQAL FHGPGALAL TEAARLFLRQ ERDARQLRVC PAAYGEPLQA AASALGAAVR LHLELGQPAA AAALCLELAA ALRDLGQPAA AAGHFQRAAQ LQLPQLPLAA LQALGEAASC QLLARDYTGA LAVFTRMQL AREHGSHPVQ SLPPPPPPAP QPGPGATPAL

image



Background

Alternative Names 40-kDa huntingtin-associated protein (HAP40) (CpG island protein) (Factor VIII intron 22 protein)

Note

Product Name: Recombinant Human F8A3,F8A1,F8A2 Protein
Catalog #: PHH10052



For Research Use Only , Not for Diagnostic Use