



Product Name: Recombinant Human PDZD9 Protein
Catalog #: PHH10148



Summary

Name	PDZD9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	Q8IXQ8
Host	HEK293
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PVTSTPKKIE LAKDESFTSS DDNENVLDLK RLQYYRYPWS TVHHPARRPI SISRDPWHGYK KKNHTISVGK DINCDVMIHR DDKKEVRAPS PYWIMVKQDN ESSSSSTSST SDAFWLEDCA QVEEGKAQLV

image



Background

Alternative Names PDZ domain-containing protein 9

Note

For Research Use Only , Not for Diagnostic Use