

Product Name: Recombinant Human ITIH3 (C-6His)
Catalog #: PHH1022




Summary

Name	ITIH3/Inter-alpha-trypsin inhibitor heavy chain H3 Protein
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Inter-alpha-trypsin Inhibitor Heavy Chain H3 is produced by our Mammalian expression system and the target gene encoding Leu35-Asp651 is expressed with a 6His tag at the C-terminus.
Accession #	Q06033
Host	Human Cells
Species	Human
Predicted Molecular Mass	70.4 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human ITIH3 (C-6His)
Catalog #: PHH1022



Alternative Names

Inter-alpha-trypsin inhibitor heavy chain H3;ITI heavy chain H3;ITI-HC3;Inter-alpha-inhibitor heavy chain 3;Serum-derived hyaluronan-associated protein;SHAP

Background

ITIH3, which is short for Inter-alpha-trypsin inhibitor heavy chain H3, is a 890 aa. protein. It is secreted expression, and belongs to the ITIH family. I-alpha-I plasma protease inhibitors are assembled from one or two heavy chains (H1, H2 or H3) and one light chain, bikunin. Inter-alpha-inhibitor (I-alpha-I) is composed of H1, H2 and bikunin, inter-alpha-like inhibitor (I-alpha-LI) of H2 and bikunin, and pre-alpha-inhibitor (P-alpha-I) of H3 and bikunin. ITIH3 may act as a carrier of hyaluronan in serum or as a binding protein between hyaluronan and other matrix protein, including those on cell surfaces in tissues to regulate the localization, synthesis and degradation of hyaluronan which are essential to cells undergoing biological processes.

Note

For Research Use Only , Not for Diagnostic Use.