



Product Name: Recombinant Human S100A6 Protein
Catalog #: PHH10496



Summary

Name	S10A6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	P06703
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLVAIFHKYS GREGDKHTLS KKEKELIQK ELTIGSKLQD AEIARLMEDL DRNKDQEVNF QEYVTFGLGAL ALIYNEALKG

image



Background

Alternative Names Protein S100-A6 (Calcyclin) (Growth factor-inducible protein 2A9) (MLN 4) (Prolactin receptor-associated protein) (PRA) (S100 calcium-binding protein A6)

Note

For Research Use Only , Not for Diagnostic Use