

Product Name: Recombinant Human CDK5R1 Protein
Catalog #: PHH10515



Summary

Name	CD5R1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15078
Host	HEK293
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YRKATLFEDG AATVGHYTAV QNSKNAKDKN LKRHSIISVL PWKRIVAVSA KKKNSKKVQP NSSYQNNITH LNNENLKKSL SCANLSTFAQ PPPAQPPAPP ASQLSGSQTG GSSSVKKAPH PAVTSAGTPK RVIVQASTSE

image



Background

Alternative Names	Cyclin-dependent kinase 5 activator 1 (CDK5 activator 1) (Cyclin-dependent kinase 5 regulatory subunit 1) (TPKII regulatory subunit) [Cleaved into: Cyclin-dependent kinase 5 activator 1, p35 (p35); Cyclin-dependent kinase 5 activator 1, p25 (p25) (Tau protein kinase II 23 kDa subunit) (p23)]
--------------------------	--

Note

Product Name: Recombinant Human CDK5R1 Protein
Catalog #: PHH10515



For Research Use Only , Not for Diagnostic Use