


Product Name: Recombinant Human MAP3K11 Protein
Catalog #: PHH10519



Summary

Name	M3K11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16584
Host	HEK293
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEELTRAARE QRSQAEQLRR REHLLAQWEL EVFERELTLL LQQVDRERPH VRRRRGTGFKR SKLRARDGGE RISMPDFFKH RITVQASPLG DRRRNVEFVG PGDSPTFPRF RAIQLEPAEP GQAWGRQSPR RLEDSSNGER RACWAWGPSS PKPGEAQNGR RRSRMDEATW YLDSDDSSPL

image



Background

Alternative Names	Mitogen-activated protein kinase kinase kinase 11 (EC 2.7.11.25) (Mixed lineage kinase 3) (Src-homology 3 domain-containing proline-rich kinase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MAP3K11 Protein
Catalog #: PHH10519

