


Product Name: Recombinant Human DYRK2 Protein
Catalog #: PHH10578




Summary

Name	DYRK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q92630
Host	HEK293
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QQLFEDNSNK RTVLTTQPNG LTTVGKTGLP VVPERQLDSI HRRQGSSTSL KSMEGMGKVK ATPMTPEQAM KQYMQKLTAH EHHEIFSYPE IYFLGLNAKK RQGMTGGPNN GGYDDDQGSY VQVPHDHVAY RYEV LKLVIGK

image



SDS



WB

Background

Alternative Names Dual specificity tyrosine-phosphorylation-regulated kinase 2 (EC 2.7.12.1)

Note

For Research Use Only , Not for Diagnostic Use