



Product Name: Recombinant Human SNRPG Protein
Catalog #: PHH10680



Summary

Name	RUXG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	P62308
Host	HEK293
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KFMDKKLSLK LNGGRHVQGI LRGFDPFMNL VIDECEMAT SGQQNNIGMV VIRGNSIIML

image



Background

Alternative Names Small nuclear ribonucleoprotein G (snRNP-G) (Sm protein G) (Sm-G) (SmG)

Note

For Research Use Only , Not for Diagnostic Use