



Product Name: Recombinant Human INA Protein
Catalog #: PHH10705



Summary

Name	AINX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q16352
Host	HEK293
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QYESLAAKNL QSAEEWYKSK FANLNEQAAR STEAIRASRE EIHEYRRQLQ ARTIEIEGLR GANESLERQI LELEERHSAE VAGYQDSIGQ LENDLRNTKS

image



Background

Alternative Names	Alpha-internexin (Alpha-Inx) (66 kDa neurofilament protein) (NF-66) (Neurofilament-66) (Neurofilament 5)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use