

Product Name: Recombinant Human Lefty-A (N-6His)
Catalog #: PHH1071




Summary

Name	Left-right Determination Factor 2
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Left-right Determination Factor 2 is produced by our Mammalian expression system and the target gene encoding Phe78-Pro366 is expressed with a 6His tag at the N-terminus.
Accession #	O00292
Host	Human Cells
Species	Human
Predicted Molecular Mass	39.1 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20 mM His-HCl, 10% Trehalose, 4% Mannitol, 100 mM NaCl, 0.1% Tween 80, pH5.5.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human Lefty-A (N-6His)
Catalog #: PHH1071



Alternative Names

Left-right determination factor 2; Endometrial bleeding-associated factor; Left-right determination factor A; Protein lefty-2; Protein lefty-A; Transforming growth factor beta-4; TGF-beta-4; LEFTY2; EBAF; LEFTA; LEFTYA; TGFβ4

Background

Left-right determination factor 2 (LEFTY2) is a secreted protein which belongs to the TGF-beta family. Lefty was first identified in a screen for undifferentiated cell-specific cDNAs from the P19 mouse embryonal carcinoma cells. Its mRNA expression on the left side of the developing embryo earned the name "Lefty". The human orthologue was initially identified as EbaF, Endometrial bleeding associated factor. Lefty contains the six cysteine residues that are conserved among TGF-β related proteins and that are necessary to form the cysteine knot structure. Its function in patterning left-right asymmetry of the developing organ systems such as the heart and lung is consistent in all vertebrate species examined. Lefty acts as an antagonist to Nodal signaling, potentially by competing for binding to a common receptor. It may play a role in endometrial bleeding.

Note

For Research Use Only , Not for Diagnostic Use.