



Product Name: Recombinant Human CCDC25 Protein
Catalog #: PHH10763



Summary

Name	CCD25_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q86WR0
Host	HEK293
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDIPKEVLMD CAHLVKANSI QGCKMNNVNV VYTPWSNLKK TADMDVGQIG FHRQKDVKIV TVEKKVNEIL NRLEKTKVER FPDLAAEKEC RDREERNEKK AQIQEMKKRE KEEMKKKREM DELRSYSSLM KVENMSSNQD

image



Background

Alternative Names Coiled-coil domain-containing protein 25

Note

For Research Use Only , Not for Diagnostic Use