



Product Name: Recombinant Human ADGRL3 Protein
Catalog #: PHH10908



Summary

Name	AGRL3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9HAR2
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPSVKDISTT GPLGMGSTTT STTLRTTTLT PGRSTTPSVS GRRNRSTSTP SPAVEVLDDM TTHLPSASSQ IPALEESCEA

image



Background

Alternative Names	Adhesion G protein-coupled receptor L3 (Calcium-independent alpha-latrotoxin receptor 3) (CIRL-3) (Latrophilin-3) (Lectomedin-3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use