


**Product Name: Recombinant Human REG3G (C-6His)**  
**Catalog #: PHH1421**



## Summary

<b>Name</b>	REG3G/Regenerating islet-derived protein 3-gamma//Reg3G/Reg III-gamma
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	<1 EU/μg as determined by LAL test.
<b>Construction</b>	Recombinant Human Regenerating Islet-derived Protein 3-gamma is produced by our Mammalian expression system and the target gene encoding Glu27-Asp175 is expressed with a 6His tag at the C-terminus.
<b>Accession #</b>	Q6UW15
<b>Host</b>	Human Cells
<b>Species</b>	Human
<b>Predicted Molecular Mass</b>	17.5 KDa
<b>Formulation</b>	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Reconstitution</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE image



## Background

**Product Name: Recombinant Human REG3G (C-6His)**  
**Catalog #: PHH1421**



---

**Alternative Names**

**Background**

Regenerating islet-derived protein 3-gamma; REG-3-gamma; Pancreatitis-associated protein 1B; REG3G; PAP-1B; Regenerating islet-derived protein III-gamma Pancreatitis-associated protein 1B, Regenerating islet-derived protein III-gamma, is a secreted protein which contains 1 C-type lectin domain. It is expressed almost exclusively in the pancreas and also expressed in testis, but not found in small intestine. This protein might be a stress protein involved in the control of bacterial proliferation.

**Note**

For Research Use Only , Not for Diagnostic Use.