

Product Name: Recombinant Human SOST (C-6His)
Catalog #: PHH1469




Summary

Name	Sclerostin/SOST
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Sclerostin is produced by our Mammalian expression system and the target gene encoding Gln24-Tyr213 is expressed with a 6His tag at the C-terminus.
Accession #	Q9BQB4
Host	Human Cells
Species	Human
Predicted Molecular Mass	22.3 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM NaAc-Hac, 8% Trehalose, 0.02% Tween 80, pH5.0.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human SOST (C-6His)
Catalog #: PHH1469



Alternative Names

Sclerostin;SOST; UNQ2976; PRO7455; PRO7476

Background

Sclerostin, also known as SOST, is a member of the Cerberus/DAN family of BMP antagonists. SOST is a secreted glycoprotein with a C-terminal cysteine knot-like (CTCK) domain. It shows sequence similarity to the DAN (differential screening-selected gene aberrative in neuroblastoma) family of bone morphogenetic protein (BMP) antagonists. SOST is produced by the osteocyte and has anti-anabolic effects on bone formation. It is a negative regulator of bone growth that acts through inhibition of Wnt signaling and bone formation.

Note

For Research Use Only , Not for Diagnostic Use.