


Product Name: Recombinant Human Serpin A8 (C-6His)
Catalog #: PHH1502



Summary

Name	Serpin A8/Angiotensinogen
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Serine Protease Inhibitor-clade A8 is produced by our Mammalian expression system and the target gene encoding Asp34-Ala485 is expressed with a 6His tag at the C-terminus.
Accession #	AAH11519.1
Host	Human Cells
Species	Human
Predicted Molecular Mass	50.8 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Product Name: Recombinant Human Serpin A8 (C-6His)
Catalog #: PHH1502



Alternative Names

Angiotensinogen; Serpin A8; AGT and SERPINA8

Background

Angiotensinogen is involved in maintaining blood pressure and in the pathogenesis of essential hypertension and preeclampsia. Genetic variations in Angiotensinogen are associated with susceptibility to essential hypertension, and can cause renal tubular dysgenesis, a severe disorder of renal tubular development. Defects in the encoding gene have also been associated with non-familial structural atrial fibrillation, and inflammatory bowel disease.

Note

For Research Use Only , Not for Diagnostic Use.