


Product Name: Recombinant Human 4-1BB (C-Fc-6His)
Catalog #: PHH1697



Summary

Name	4-1BB/CD137/TNFRSF9/ILA
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human 4-1BB Ligand Receptor is produced by our Mammalian expression system and the target gene encoding Leu24-Gln186 is expressed with a human IgG1 Fc, 6His tag at the C-terminus.
Accession #	Q07011
Host	Human Cells
Species	Human
Predicted Molecular Mass	44 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Product Name: Recombinant Human 4-1BB (C-Fc-6His)
Catalog #: PHH1697



Alternative Names

CD137; ILA; TNFRSF9; 4-1BB ligand receptor; CDw137; T-cell antigen 4-1BB homolog; T-cell antigen ILA

Background

Tumor necrosis factor receptor superfamily member 9 (TNFRSF9) is an inducible T cell surface protein belonging to the TNF receptor superfamily. It is a single-pass type I membrane protein which contains 4 TNFR-Cys repeats. The human and mouse proteins share 60% amino acid sequence identity. It is absent from naive T cells, but upregulated and continually expressed following T cell activation. It is a receptor for TNFSF9/4-1BBL, and possibly active during T cell activation.

Note

For Research Use Only , Not for Diagnostic Use.