


Summary

| | |
|---------------------------------|--|
| Name | TREML1/TLT-1/Trem-like transcript 1 protein |
| Purity | Greater than 95% as determined by reducing SDS-PAGE |
| Endotoxin level | <1 EU/μg as determined by LAL test. |
| Construction | Recombinant Human Triggering Receptor Expressed on Myeloid Cells-Like Protein 1 is produced by our Mammalian expression system and the target gene encoding Gln16-Pro162 is expressed with a human IgG1 Fc tag at the C-terminus. |
| Accession # | Q86YW5 |
| Host | Human Cells |
| Species | Human |
| Predicted Molecular Mass | 42.9 KDa |
| Formulation | Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4. |
| Shipping | The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Reconstitution | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |

SDS-PAGE image



Product Name: Recombinant Human TREML1 (C-Fc)
Catalog #: PHH1720



Background

Alternative Names

Trem-Like Transcript 1 Protein; TLT-1; Triggering Receptor Expressed on Myeloid Cells-Like Protein 1; TREML1; TLT1

Background

Triggering Receptor Expressed on Myeloid Cells-Like Protein 1 (TREML1) is a single-pass type I membrane protein. TREML1 precursor contains a 15 amino acid signal peptide, a 147 amino acid extracellular domain with an Ig-like V-type (immunoglobulin-like) domain, and 128 amino acid cytoplasmic domain. It can be expressed exclusively in platelets and megakaryocytes (MKs). It is a cell surface receptor that may play a role in the innate and adaptive immune response. TREML1 is sequestered in cytoplasmic vesicles in resting platelets. TREML1 is transported to the cell surface after stimulation by thrombin. Soluble fragments can be released into the serum by proteolysis. The phosphorylated TREML1 can interact with PTPN6 and PTPN11. TREML1 may participate in maintaining vascular hemostasis and regulating coagulation and inflammation at sites of injury.

Note

For Research Use Only , Not for Diagnostic Use.