

Product Name: Recombinant Human FAM172A (C-6His)
Catalog #: PHH1864




Summary

Name	FAM172A
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Protein FAM172A/FAM172A is produced by our Mammalian expression system and the target gene encoding Gln19-Leu416 is expressed with a 6His tag at the C-terminus.
Accession #	Q8WUF8
Host	Human Cells
Species	Human
Predicted Molecular Mass	47 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human FAM172A (C-6His)
Catalog #: PHH1864



Alternative Names

Protein FAM172A;C5orf21

Background

FAM172A is encoded by FAM172A gene. It is a 416 aa. protein, and belongs to the UPF0528 family. This protein has 4 isoforms produced by alternative splicing. There is a natural variant location which is S131N. FAM172A has a 18 aa. signal peptide which helps it to be secreted protein.

Note

For Research Use Only , Not for Diagnostic Use.