


Summary

Name	Fc gamma RI/CD64/FCGR1A/Fc γ RI
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/ μ g as determined by LAL test.
Construction	Recombinant Human High Affinity Immunoglobulin Gamma Fc Receptor I is produced by our Mammalian expression system and the target gene encoding Gln16-Pro288 is expressed with a 6His tag at the C-terminus.
Accession #	P12314
Host	Human Cells
Species	Human
Predicted Molecular Mass	31.7 KDa
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8 $^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μ g/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Alternative Names

Product Name: Recombinant Human CD64 (C-6His)
Catalog #: PHH2073



High affinity immunoglobulin gamma Fc receptor I; IgG Fc receptor I; Fc-gamma RI; FcRI; Fc-gamma RIA; FcgammaRIa; CD64; FCGR1A

Background

CD64 (FcγRI), one of the Fc receptors for IgG, is a membrane glycoprotein that mediates endocytosis, phagocytosis, antibody-dependent cellular cytotoxicity, cytokine release, and superoxide production. CD64 is also structurally distinct, containing an extracellular Ig-interactive region of three Ig-like domains in contrast to the two domains of the low affinity receptors FcγRII and FcγRIII. It is normally expressed on the surfaces of monocytes and macrophages.

Note

For Research Use Only , Not for Diagnostic Use.