

Product Name: Recombinant Human B4GALT1 (C-6His)
Catalog #: PHH2132




Summary

Name	B4GALT1
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Beta-1,4-galactosyltransferase 1 is produced by our Mammalian expression system and the target gene encoding Gly44-Ser398(Tyr285Leu) is expressed with a 6His tag at the C-terminus.
Accession #	P15291
Host	Human Cells
Species	Human
Predicted Molecular Mass	40.1 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 2mM EDTA, 20% Glycerol, pH 8.0.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image

Product Name: Recombinant Human B4GALT1 (C-6His)
Catalog #: PHH2132



Alternative Names

1,4- galactosyltransferase, polypeptide 1; B4GalT1; B4GAL-T1; beta-1,4-galactosyltransferase 1; Beta-1,4-GalTase 1; beta4Gal-T1; betaGlcNAc beta; CDG2D; GT1; GTB

Background

Beta1,4-Galactosyltransferase-I (B4GALT1), one of seven beta1,4-galactosyltransferases, is an enzyme commonly found in the trans-Golgi complex that adds galactose to oligosaccharides. By sequence similarity, the beta 4GalTs form four groups: beta 4GalT1 and beta 4GalT2, beta 4GalT3 and beta 4GalT4, beta 4GalT5 and beta 4GalT6, and beta 4GalT7. beta 4GalT1 is unique among the seven enzymes because it can be expressed either as membrane associated form or secreted form. The secreted form is restricted to lactating mammary tissues where the enzyme forms a heterodimer with alpha -lactalbumin to catalyze the synthesis of lactose. The Golgi complex form catalyzes the production of lactose in the lactating mammary gland and could also be responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids. The cell surface form functions as a recognition molecule during a variety of cell to cell and cell to matrix interactions, as those occurring during development and egg fertilization, by binding to specific oligosaccharide ligands on opposing cells or in the extracellular matrix. Defects in beta 4GalT1 are the cause of congenital disorder of glycosylation type 2D (CDG2D) .

Note

For Research Use Only , Not for Diagnostic Use.