


Product Name: Recombinant Mouse BTNL9 (C-Fc)
Catalog #: PHH2309



Summary

Name	BTNL9/Butyrophilin-like Protein 9
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Mouse Butyrophilin-like Protein 9 is produced by our Mammalian expression system and the target gene encoding Asp36-Lys257 is expressed with a human IgG1 Fc tag at the C-terminus.
Accession #	Q8BJE2
Host	Human Cells
Species	Human
Predicted Molecular Mass	51.8 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Background

Product Name: Recombinant Mouse BTNL9 (C-Fc)
Catalog #: PHH2309



Alternative Names

Butyrophilin-like protein 9; BTNL9

Background

Butyrophilin-Like Protein 9 (BTNL9) is single-pass type I membrane protein member of the BTN/MOG family that belongs to the immunoglobulin superfamily. BTNL9 consists of two domains: one B30.2/SPRY domain and one Ig-like V-type (immunoglobulin-like) domain. Human BTNL9 mRNA has been identified in adipose, lung, thymus, spleen, colon, and cardiac tissues, but its highest levels of expression were found in B cells. BTNL9 expression has also been found to be down-regulated in colon cancer tumors.

Note

For Research Use Only , Not for Diagnostic Use.