


Product Name: Recombinant Human TNFRSF11B Protein
Catalog #: PHH2705



Summary

Name	TR11B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O00300
Host	HEK293
Species	Human
Molecular Mass	86 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ETFPPKYLHY DEETSHQLLC DKCPPGTLYK QHCTAKWKTVCAPCPDHYYT DSWHTSDECL YCSPVCKELQ YVKQECNRTH NRVCECKEGR YLEIEFCLKH RSCPPGFGVV QAGTPERNTV CKRCPDGFFS NETSSKAPCR KHTNCSVFGL LLTQKGNATH DNICSGNSES TQKCGIDVTL CEEAFFRFAV PTKFTPNWLS VLVDNLPGTK VNAESVERIK RQHSSQEQTF QLLKLWKHQNKDQDIVKKII QDIDLCENSV QRHIGHANLT FEQLRSLMES LPGKKVGAED IEKTIKACKP SDQILKLLSL WRIKNGDQDT LKGLMHALKH SKTYHFPKTV TQSLKKTIRF LHSFTMYKLY QKLFLEMIGN QVQSVKISCL

image



Background

Alternative Names	Tumor	necrosis	factor	receptor	superfamily	member	11B
--------------------------	-------	----------	--------	----------	-------------	--------	-----



Product Name: Recombinant Human TNFRSF11B Protein
Catalog #: PHH2705

(Osteoclastogenesis inhibitory factor) (Osteoprotegerin)

Note

For Research Use Only , Not for Diagnostic Use