


Product Name: Recombinant Human MUSK Protein
Catalog #: PHH2712



Summary

Name	MUSK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15146
Host	HEK293
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TTPLETVDAL VEEVATFMCA VESYPQPEIS WTRNKILIKL FDTRYSIREN GQLLTILSVE DSDDGIYCCT ANNGVGGAVE SCGALQVKMK PKITRPPINV KIIEGLKAVL PCTTMGNPKP SVSWIKGDSP LRENSRIAVL ESGSLRIHNV QKEDAGQYRC VAKNSLGTAY SKVVKLEVEV FARILRAPES HNVTFGSFVT LHCTATGIPV PTITWIENG N AVSSGSIQES VKDRVIDSRL QLFITKPGLY TCIATNKHGE KFSTAKAAAT ISIAEWSKPQ KDNKGYCAQY RGEVCNAVLA KDALVFLNTS YADPEEAQEL LVHTAWNELK VVSPVCRPAA EALLCNHIFQ ECSPGVVPTP IPICREYCLA VKELFCAKEW LVMEETHRG LYRSEMHLLS VPECSKLPSM HWDPTACARL PHLDYNKENL KTFPPMTSSK PSVDIPNLPS SSSSFSVSP

image



Background

Product Name: Recombinant Human MUSK Protein
Catalog #: PHH2712



Alternative Names

Muscle, skeletal receptor tyrosine-protein kinase (EC 2.7.10.1) (Muscle-specific tyrosine-protein kinase receptor) (MuSK) (Muscle-specific kinase receptor)

Note

For Research Use Only , Not for Diagnostic Use