


Product Name: Recombinant Human SLC26A4 Protein
Catalog #: PHH2725



Summary

Name	S26A4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O43511
Host	HEK293
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGSIPSTDIY KSTKNYKNIE EPQGVKILRF SSPIFYGNVD GFKKCIKSTV GFDAIRVYNK RLKALRKIQK LIKSGQLRAT KNGIISDAVS TNNAFEPDED IEDLEELDIP TKEIEIQVDW NSELPVKVNV PKVPIHSLVL DCGAISFLDV

image



Background

Alternative Names	Pendrin (Sodium-independent chloride/iodide transporter) (Solute carrier family 26 member 4)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use