



Product Name: Recombinant Human LRP4 Protein
Catalog #: PHH2740



Summary

| | |
|------------------------------|--|
| Name | LRP4_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-hFc & C-His |
| Accession # | O75096 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 50 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | CSPLDFHCDN GKCIRRSWVC DGDNDCEDDS DEQDCPPREC EEDEFPCQNG YCIRSLWHCD GDNDCGDNSD EQCDMRKCSD KEFRCSDGSC IAEHWYCDGD TDCKDGSDEE NCPSAVPAPP CNLEEFQCAY GRCILDIYHC DGDDDCGDWS DESDCSSHQP CRSGEFMCDS GLCINAGWRC DGDADCDDQS DERNCTTSMC |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Low-density lipoprotein receptor-related protein 4 (LRP-4) (Multiple epidermal growth factor-like domains 7) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LRP4 Protein
Catalog #: PHH2740

