


Summary

Name	CFAB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P00751
Host	HEK293
Species	Human
Molecular Mass	84 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TPWSLARPQG SCSLEGVEIK GGSFRLQEG QALEYVCPSPG FYPYPVQTRT CRSTGSWSTL KTQDQKTVRK AECRAIHCPR PHDFENGEYW PRSPYYNVSD EISFHCYDGY TLRGSANRTC QVNGRWSGQT AICDNGAGYC SNPGIPIGTR KVGSIYRLED SVTYHCSRGL TLRGSQRRTC QEGGSWSGTE PSCQDSFMYD TPQEVAEAFI SSLTETIEGV DAEDGHGPGE QKQRKIVLDP SGSMNIYLV DGSDSIGASN FTGAKKCLVN LIEKVASYGV KPRYGLVITYA TYPKIWVKVS EADSSNADWV TKQLNEINYE DHKLKSGTNT KKALQAVYSM MSWPDDVPPE GWNRTRHVII LMTDGLHNMG GDPITVIDEI RDLLYIGKDR KNPREDYLDV YVFGVGPLVN QVNINALASK KDNEQHVFKV KDMENLEDVF YQMIDESQSL SLCGMVWEHR KGTDYHKQPW QAKISVIRPS KGHESCMGAV VSEYFVLTA HCFTVDDKEH SIKVSVGGEK RDLEIEVVLV HPNYNINGKK EAGIPEFYDY DVALIKLKNK LKYGQTIRPI CLPCTEGTTR ALRLPPTTTC QQQKEELLPA QDIKALFVSE EEKLRTRKEV YIKNGDKKGS CERDAQYAPG YDKVKDISEV VTPRFLCTGG VSPYADPNTC RGDSGGPLIV HKRSRFIQVG VISWGVVDVC KNQKRQKQVP AHARDFHINL FQVLPWLKEK LQDEDLGLF

image

Product Name: Recombinant Human CFB Protein
Catalog #: PHH2764



Background

Alternative Names

Complement factor B (EC 3.4.21.47) (C3/C5 convertase) (Glycine-rich beta glycoprotein) (GBG) (PBF2) (Properdin factor B) [Cleaved into: Complement factor B Ba fragment; Complement factor B Bb fragment]

Note

For Research Use Only , Not for Diagnostic Use