

Product Name: Recombinant Human MBP Protein
Catalog #: PHH2793



Summary

Name	MBP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P02686
Host	HEK293
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GFLPRHRDTG ILDSIGRFFG GDRGAPKRGs GKDSHHPART AHYGSLPQKS HGRTQDENPV VHFFKNIVTP RTPPSQGKG RGLSLSRFSW GAEGQRPFG YGGRASDYKS AHKGFKGVDA QGTLSKIFKL GGRDSRSGSP

image



Background

Alternative Names	Myelin basic protein (MBP) (Myelin A1 protein) (Myelin membrane encephalitogenic protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use