



**Product Name: Recombinant Human CSF2 protein**  
**Catalog #: PHH2815**



## Summary

<b>Name</b>	CSF2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-His & C-GST
<b>Accession #</b>	P04141
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	41 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RSPSPSTQPW EHVNAIQEAR RLLNLSRDTA AEMNETVEVI SEMFDLQEPT CLQTRLELYK QGLRGSCLKL KGPLTMMASH YKQHCPTPE

## image



## Background

<b>Alternative Names</b>	Granulocyte-macrophage colony-stimulating factor (GM-CSF) (Colony-stimulating factor) (CSF) (Molgramostin) (Sargramostim)
--------------------------	---

## Note

For Research Use Only , Not for Diagnostic Use