


Summary

Name	CEAM5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P06731
Host	HEK293
Species	Human
Molecular Mass	113 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLTIESTPFN VAEGKEVLLL VHNLPQHFLG YSWYKGERVD GNRQIIGYVI GTQQATPGPA YSGREIYYPN ASLLIQNIIQ NDTGFYTLHV IKSDLVNEEA TGQFRVYPEL PKPSSSNNS KPVEDKDAVA FTCEPETQDA TYLWWWVNNQS LPVSPRLQLS NGNRTLTLFN VTRNDTASYK CETQNPVSAR RSDSVILNVL YGPDAPTISP LNTSYRSGEN LNLSCHAASN PPAQYSWFVN GTFQQSTQEL FIPNITVNNS GSYTCQAHNS DTGLNRTT VT TITVYAEPK PFITSNNSNP VEDEDAVALT CEPEIQNTTY LWWWVNNQSLP VSPRLQLSND NRTLTLTLLSVT RNDVGPYECG IQNELSVDHS DPVILNVLYG PDDPTISPSY TYRPGVNLS LSCHAASNPP AQYSWLIDGN IQQHTQELFI SNITEKNSGL YTCQANNSAS GHSRTTVKTI TVSAELPKPS ISSNNSKPVE DKDAVAFTCE PEAQNTTYLW WVNGQSLPVS PRLQLSNGNR TLTLFNVTRN DARAYVCGIQ NSVSANRSDP VTLDVLYGPD TPIISPPDSS YLSGANLNLS CHSASNPSQ YSWRINGIPQ QHTQVLFIK ITPNNGTYA CFVSNLATGR NNSIVKSITV SASGTSPGLS A

image



Product Name: Recombinant Human CEACAM5 Protein
Catalog #: PHH2847



Background

Alternative Names

Carcinoembryonic antigen-related cell adhesion molecule 5
(Carcinoembryonic antigen) (CEA) (Meconium antigen 100) (CD antigen
CD66e)

Note

For Research Use Only , Not for Diagnostic Use