


Product Name: Recombinant Human CENPB Protein
Catalog #: PHH2857



Summary

Name	CENPB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P07199
Host	HEK293
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LHFLEGGEDS DSDSEEDDE EEDDEDEDD DDEEDGDEVP VPSFGAMAY FAMVKRYLTS FPIDDRVQSH ILHLEHDLVH VTRKNHARQA

image



Background

Alternative Names Major centromere autoantigen B (Centromere protein B) (CENP-B)

Note

For Research Use Only , Not for Diagnostic Use