



Product Name: Recombinant Human DBH Protein
Catalog #: PHH2893



Summary

Name	DOPO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His & N-MBP
Accession #	P09172
Host	HEK293
Species	Human
Molecular Mass	68 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APRESLPYH IPLDPEGSLE LSWNVSYTQE AIHFQLLVRR LKAGVLFGMS DRGELENADL VVLWTDGDTA YFADAWSDQK GQIHLDPQQD YQLLQVQRTP EGLTLLFKRP FGTCDPKDYI IEDGTVHLVY GILEEPPFRSL EAINGSGLQM GLQRVQLLKP NIPEPELPSD ACTMEVQAPN IQIPSQETTY WCYIKELPKG FSRHIIKYE PIVTKGNEAL VHHMEVFQCA

image



Background

Alternative Names Dopamine beta-hydroxylase (EC 1.14.17.1) (Dopamine beta-monooxygenase) [Cleaved into: Soluble dopamine beta-hydroxylase]

Note

Product Name: Recombinant Human DBH Protein
Catalog #: PHH2893



For Research Use Only , Not for Diagnostic Use