


Product Name: Recombinant Human RO60 Protein
Catalog #: PHH2911



Summary

Name	RO60_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P10155
Host	HEK293
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MNQRVLGSIL NASTVAAAMC MVTTRTEKDS YVVAFSDEMV PCPVTDDMTL QQVLMAMSQI PAGGTDCSLP MIWAQKTNTN ADVFIVFTDN ETFAGGVHPA IALREYRKKM DIPAKLIVCG MTSNGFTIAD PDDRGMLDMC GFDTGALDVI

image



Background

Alternative Names	60 kDa SS-A/Ro ribonucleoprotein (60 kDa Ro protein) (60 kDa ribonucleoprotein Ro) (Ro60) (RoRNP) (Ro 60 kDa autoantigen) (Ro60 autoantigen) (Sjogren syndrome antigen A2) (Sjogren syndrome type A antigen) (SS-A) (TROVE domain family member 2)
--------------------------	--

Note

Product Name: Recombinant Human RO60 Protein
Catalog #: PHH2911



For Research Use Only , Not for Diagnostic Use