


Product Name: Recombinant Human FGFR1 Protein
Catalog #: PHH2921



Summary

Name	FGFR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11362
Host	HEK293
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DDDSSEEKE TDNTPNRMP VAPYWTSPEK MEKKLHAVPA AKTVKFKCPS SGTPNPTRLRW LKNGKEFKPD HRIGGYKVRY ATWSIIMDSV VPSDKGNYTC IVENEYGSIN

image



SDS



WB

Background

Alternative Names	Fibroblast growth factor receptor 1 (FGFR-1) (EC 2.7.10.1) (Basic fibroblast growth factor receptor 1) (BFGFR) (bFGF-R-1) (Fms-like tyrosine kinase 2) (FLT-2) (N-sam) (Proto-oncogene c-Fgr) (CD antigen CD331)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FGFR1 Protein
Catalog #: PHH2921

