



Product Name: Recombinant Human IL10 Protein
Catalog #: PHH2982



Summary

Name	IL10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P22301
Host	HEK293
Species	Human
Molecular Mass	61 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPGQGTQSEN SCTHFPGNLP NMLRDLRDAF SRVKTFQMK DQLDNLLLKE SLLEDFKGYL GCQALSEMIQ FYLEEVMPPQA ENQDPDIKAH VNSLGENLKT LRLRLRRCHR FLPCENKSKA VEQVKNAFNK LQEKGIYKAM SEFDIFINYI EAYMTMKIRN

image



Background

Alternative Names Interleukin-10 (IL-10) (Cytokine synthesis inhibitory factor) (CSIF)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IL10 Protein
Catalog #: PHH2982

