



Summary

Name	CO4A5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P29400
Host	HEK293
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGIKGEKGER GFPGLEHHPG LPGFPGPEGP PGPRGQKGDD GIPGPPGPKG IRGPPGLPGF PGTPGLPGMP GHDGAPGPQG IPGCNGTKGE RGFPGSPGFP GLQGPPGPPG IPGMKGEPGS IIMSSLPGPK GNPGYGPPG IQGLPGPTGI PGPIGPPGPP GLMGPPGPPG LPGPKGNMGL NFQGPKGEKG EQGLQGPPGP PGQISEQKRP IDVEFQKGDQ GLPGDRGPPG PPGIRGPPGP PGGEKGEKGE QGEPGKRGKP GKDGENGQPG IPGLPGDPGY PGEPGRDGEK GQKGDTPPPG PPGLVIPRPG TGITIGEKN IGLPGLPGEK GERGFPIQG PPGLPGGPA AVMGPPGPPG FPGERGQKGD EGPPGISIPG PPGLDGQPGA PGLPGPPGA GPHIPPSDEI CEPGPPGPPG SPGDKGLQGE QGVKGDKGDT CFNCIGTGIS GPPGQPLPG LPGPPSLGF PGQKGEKGQA GATGPKGLPG IPGAPGAPG PGSKGEPGDI

image



Background

Product Name: Recombinant Human COL4A5 Protein
Catalog #: PHH2997



Alternative Names

Collagen alpha-5(IV) chain

Note

For Research Use Only , Not for Diagnostic Use