


Summary

Name	SPB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P29508
Host	HEK293
Species	Human
Molecular Mass	87 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MNSLSEANTK FMFDLFQQFR KSKENNIFYS PISITSALGM VLLGAKDNTA QQIKKVLHFD QVTENTTGKA ATYHVDRSGN VHHQFQKLLT EFNKSTDAYE LKIANKLFG E KTYLFLQEYL DAIKKFYQTS VESVDFANAP EESRKKINSW VESQTNEKIK NLIPEGNIGS NTTLVLVNAI YFKGQWEKKF NKEDTKEEF WPNKNTYKSI QMMRQYTSFH FASLEDVQAK VLEIPYKGKD LSMIVLLPNE IDGLQKLEEK LTAEKLMWT SLQNMRETRV DLHLPRFKVE ESYDLKDTLR TMGMVDIFNG DADLSGMTGS RGLVLSGVLH KAFVEVTEEG AAAAAATAVV GFGSSPTSTN EEHCNHPFL FFIRQNKTNS ILFYGRFSSP

image



Background

Alternative Names Serpin B3 (Protein T4-A) (Squamous cell carcinoma antigen 1) (SCCA-1)

Product Name: Recombinant Human SERPINB3 Protein
Catalog #: PHH3000



Note

For Research Use Only , Not for Diagnostic Use