


Product Name: Recombinant Human PFN2 Protein
Catalog #: PHH3016



Summary

Name	PROF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P35080
Host	HEK293
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGWQSYVDNL MCDGCCQEAA IVGYCDAYV WAATAGGVFQ SITPIEDMI VGKDREGFFT NGLTLGAKKC SVIRDSLYVD GDCTMDIRTK SQGGPEYTNV AVGRAGRVLV FVMGKEGVHG GGLNKKAYSM AKYLRDSGF

image



SDS



WB

Background

Alternative Names Profilin-2 (Profilin II)

Note

For Research Use Only , Not for Diagnostic Use