


## Summary

<b>Name</b>	TPO_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-His & N-SUMO
<b>Accession #</b>	P40225
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	47 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SPAPPACDLR VLSKLLRDSH VLHSRLSQCP EVHPLPTPVL LPAVDFSLGE WKTQMEETKA QDILGAVTLL LEGVMAARGQ LGPTCLSSLL GQLSGQVRLL LGALQSLLGT QLPPQGRRTA HKDPNAIFLS FQHLLRGKVR FLMLVGGSTL CVRRAPPTTA VPSRTSLVLT LNELPNRTSG LLETNFTASA RTTGSGLLKW QQGFRAKIPG LLNQTSLRSLD QIPGYLNRIH ELLNGTRGLF PGPSRRTLGA PDISSGTSDT GSLPPNLQPG YSPSPHPPT GQYTLFPLPP TLPTPVVQLH PLLPDPSAPT PTPTSPLLNT SYTHSQNLSQ EG

## image



## Background

<b>Alternative Names</b>	Thrombopoietin (C-mpl ligand) (ML) (Megakaryocyte colony-stimulating factor) (Megakaryocyte growth and development factor) (MGDF)
--------------------------	---

**Product Name: Recombinant Human THPO Protein**  
**Catalog #: PHH3040**



---

(Myeloproliferative leukemia virus oncogene ligand)

**Note**

For Research Use Only , Not for Diagnostic Use