


Product Name: Recombinant Human GRIA1 Protein
Catalog #: PHH3049



Summary

Name	GRIA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P42261
Host	HEK293
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AATDAQAGGD NSSVQRNTYI VTTILEDPYV MLKKNANQFE GNDRYEGYCV ELAAEIAKHV GYSYRLEIVS DGKYGARDPD TKAWNGMVGE

image



Background

Alternative Names Glutamate receptor 1 (GluR-1) (AMPA-selective glutamate receptor 1) (GluR-A) (GluR-K1) (Glutamate receptor ionotropic, AMPA 1) (GluA1)

Note

For Research Use Only , Not for Diagnostic Use