



**Product Name: Recombinant Human GRIA3 Protein**  
**Catalog #: PHH3052**



## Summary

<b>Name</b>	GRIA3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-hFc & C-His
<b>Accession #</b>	P42263
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	NEYIEQRKPC DTMKVGGNLD SKGYGVATPK GSALRNAVNL AVLKLNEQGL LDKLNKWWY DKGECGSGGG

## image



## Background

**Alternative Names** Glutamate receptor 3 (GluR-3) (AMPA-selective glutamate receptor 3) (GluR-C) (GluR-K3) (Glutamate receptor ionotropic, AMPA 3) (GluA3)

## Note

For Research Use Only , Not for Diagnostic Use