



Summary

Name	COMP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49747
Host	HEK293
Species	Human
Molecular Mass	82 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGQSPLGSDL GPQMLRELQE TNAALQDVRE LLRQQVREIT FLKNTVMECD ACGMQQSVRT GLPSVRPLLH CAPGFCFPGV ACIQTESGAR CGPCPAGFTG NGSHCTDVNE CNAHPCFPRV RCINTSPGFR CEACPPGYSG PTHQGVGLAF AKANKQVCTD INECETGQHN CVPNSVCINT RGSFQCGPCQ PGFVGDQASG CQRRARFCP DGSPSECHEH ADCVLERDGS RSCVCAVGWA GNGILCGRDT DLDGFPDEKL RCPERQCRKD NCVTVPNSGQ EDVDRDGIGD ACDPADGDG VPNEKDNCPV VRNPDQRNTD EDKWGDACDN CRSQKNDDQK DTDQDGRGDA CDDDIDGDRI RNQADNCPRV PNSDQKSDG DGIGDACDNC PQKSNPDQAD VDHDFVGDAC DSDQDQDGDG HQDSRDNCPT VPNSAQEDSD HDGQGDACDD DDDNDGVPDS RDNCRLVPPN GQEDADRDRGV GDVCQDDFDA DKVVDKIDVC PENAEVTLTD FRAFQTVVLD PEGDAQIDPN WVLNQGREGI VQTMNSDPGL AVGYTAFNGV DFEGTFHVNT VTDDDYAGFI FGYQDSSSFY VVMWKQMEQT YWQANPFRV AEPGIQLKAV KSSTGPGEQL RNALWHTGDT ESQVRLWLD PRNVGWKDKK SYRWFLQHRP QVGYIRVRFY EGPELVADSN VVLDTTMRGG RLGVFCFSQE NIIWANLRYR CNDTIPEDYE THQLRQA

image

Product Name: Recombinant Human COMP Protein
Catalog #: PHH3068



Background

Alternative Names

Cartilage oligomeric matrix protein (COMP) (Thrombospondin-5) (TSP5)

Note

For Research Use Only , Not for Diagnostic Use