



Product Name: Recombinant Human CDH4 Protein
Catalog #: PHH3082



Summary

Name	CADH4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P55283
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VYNGSVDEGS KPGTYVMTVT ANDADDSTTA NGMVRYRIVT QTPQSPSQNM FTINSETGDI VTVAAGLDRE KVQQYTVIVQ ATDMEGNLNY GLSNTATAII TVTDVNDNPP

image



Background

Alternative Names Cadherin-4 (Retinal cadherin) (R-CAD) (R-cadherin)

Note

For Research Use Only , Not for Diagnostic Use